




Attaching ODFS[®] Pace to the Leg Strap

Joe Green, Senior Physiotherapist
Laura Humphreys, Physiotherapy Lead

The ODFS[®] Pace case comes with a neat system to enable attachment to the leg strap.



To remove the belt clip from the pouch, use a biro to press down onto the release button on the pouch and slide the clip in the direction of the yellow arrow above.



To attach the pouch to the leg strap, simply slide the connections together in the direction of the green arrows. To remove the pouch simply press down the button identified above with the yellow arrow, whilst sliding the pouch upwards

Where a user requires an even closer fit to the leg, the strap can be fed through the belt loop on the pouch, as in the pictures below.




Donning and Doffing a Leg Strap

Where fine hand dexterity is a problem, patients should avoid un-threading the leg strap fully when removing it, simply loosening it to go over the foot. This strategy will mean that the end of the strap does not have to be threaded through the belt loop or buckle, on a daily basis. In the photo below the strap has been loosened only to allow it to be passed over the foot.



Attaching a Footswitch without Additional Inner Sole.

In some instances, it is not appropriate or desirable to insert an additional inner sole into footwear (sandals, snug fitting shoes). In these cases the use of Velcro dots can be useful to allow firm attachment of the switch directly to the footwear, while maintaining the ability to move the switch from shoe to shoe.



Velcro Dot (hooks) on flat surface of footswitch

Velcro Dot (loops)